

Product datasheet

Anti-LI Cadherin antibody ab37788

1 图像

概述

产品名称	Anti-LI Cadherin抗体
描述	鸡多克隆抗体to LI Cadherin
宿主	Chicken
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Synthetic peptide, corresponding to N terminal amino acids 23 - 146 of Human LI Cadherin

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Constituent: PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgY

应用

Our [Abpromise guarantee](#) covers the use of **ab37788** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/2000. Predicted molecular weight: 92 kDa.

靶标

**功能** Cadherins are calcium dependent cell adhesion proteins. They preferentially interact with themselves in a homophilic manner in connecting cells; cadherins may thus contribute to the

sorting of heterogeneous cell types. LI-cadherin may have a role in the morphological organization of liver and intestine. Involved in intestinal peptide transport.

**组织特异性**

Expressed in the gastrointestinal tract and pancreatic duct. Not detected in kidney, lung, liver, brain, adrenal gland and skin.


**序列相似性**

Contains 7 cadherin domains.

**细胞定位**

Cell membrane.

**图片**



Anti-LI Cadherin antibody (ab37788) at 1/2000 dilution + E. coli derived fusion protein at 2.5 µg

**Secondary**

HRP conjugated Goat anti IgY at 1/1000 dilution

**Predicted band size:** 92 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

**Our Abpromise to you: Quality guaranteed and expert technical support**

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

**Terms and conditions**

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors