

Product datasheet

Anti-LRRC14B antibody ab122435

2 图像

概述

产品名称	Anti-LRRC14B抗体
描述	兔多克隆抗体to LRRC14B
宿主	Rabbit
经测试应用	适用于: IHC-P, ICC/IF
种属反应性	与反应: Human
免疫原	antigen sequence: LRCIEFPVPK DCYPEGAAYP QDELAMSKFN QQKYDEIAEE LRAVLLRADR EDIQVSTPLF GSFDPDIQET SNELGAFLLQ AFKTAL, corresponding to C terminal amino acids 418-503 of Human LRRC14B protein (UniProt ID: AGNHZ5). <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>
阳性对照	Human bronchus tissues

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 59% PBS, 40% Glycerol
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab122435** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
IHC-P		1/50 - 1/200. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.


应用	Ab评论	说明
----	------	----

ICC/IF      Use a concentration of 1 - 4 µg/ml.  
 Recommend PFA Fixation and Triton X-100 treatment

**靶标**


**序列相似性**      Belongs to the PRAME family.  
 Contains 5 LRR (leucine-rich) repeats.

**图片**



Immunofluorescent staining of Human cell line A-431 shows positivity in nucleoli. Recommended concentration of ab122435 1-4 µg/ml. Cells treated with PFA/Triton X-100.

Immunocytochemistry/ Immunofluorescence - Anti-LRRC14B antibody (ab122435)



ab122435, at a 1/50 dilution, staining LRRC14B in paraffin embedded Human bronchus tissue by Immunohistochemistry.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-LRRC14B antibody (ab122435)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors