

## Product datasheet

# Anti-LYRM1 antibody ab106665

### 3 图像

#### 概述

产品名称	Anti-LYRM1抗体
描述	兔多克隆抗体to LYRM1
宿主	Rabbit
特异性	Two isoforms of LYRM1 are known to exist; ab106665 will recognize both isoforms. ab106665 is predicted not to cross-react with other LYRM family members.
经测试应用	适用于: IHC-P, ICC/IF, WB
种属反应性	与反应: Mouse, Rat, Human
免疫原	A 17 amino acid synthetic peptide near the center of Human LYRM1 (EAW66844).
阳性对照	Human liver tissue lysate; Rat liver tissue.

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at 4°C (stable for up to 12 months).
存储溶液	Preservative: 0.02% Sodium azide Constituent: PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

#### 应用

Our [Abpromise guarantee](#) covers the use of **ab106665** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明
IHC-P		Use a concentration of 5 µg/ml.
ICC/IF		Use a concentration of 20 µg/ml.

应用	Ab评论	说明
WB		Use a concentration of 1 - 2 µg/ml. Predicted molecular weight: 28 kDa.

## 靶标


功能	May promote cell proliferation and inhibition of apoptosis of preadipocytes.
组织特异性	High levels in adipose tissue.
序列相似性	Belongs to the complex I LYR family.
细胞定位	Nucleus.

## 图片



Immunofluorescence of LYRM1 in rat liver tissue with ab106665 at 20 µg/mL.

Immunocytochemistry/ Immunofluorescence - Anti-LYRM1 antibody (ab106665)



**Lane 1** : Anti-LYRM1 antibody (ab106665) at 1 µg/ml


**Lane 2** : Anti-LYRM1 antibody (ab106665) at 2 µg/ml

**All lanes** : Human liver tissue lysate

Lysates/proteins at 15 µg per lane.

**Predicted band size:** 28 kDa

Western blot - Anti-LYRM1 antibody (ab106665)



ab106665, at 5µg/ml, staining LYRM1 in rat liver tissue by Immunohistochemistry.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-LYRM1 antibody (ab106665)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors