

Product datasheet

Anti-Leukotriene A4 hydrolase antibody [EPR5712] ab109434

重组 RabMAb

3 图像

概述


产品名称	Anti-Leukotriene A4 hydrolase 抗体[EPR5712]
描述	兔单克隆抗体[EPR5712] to Leukotriene A4 hydrolase
宿主	Rabbit
经测试应用	适用于: WB, IP, IHC-P, Flow Cyt 不适用于: ICC
种属反应性	与反应: Rat, Human
免疫原	Synthetic peptide within Human Leukotriene A4 hydrolase aa 1-100 (internal sequence). The exact sequence is proprietary.
阳性对照	HeLa, 293T, T47-D, and A549 cell lysates; Human breast carcinoma tissue
常规说明	Mouse: We have preliminary internal testing data to indicate this antibody may not react with this species. Please contact us for more information.

Our RabMAb[®] technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to [RabMAb[®] patents](#)

This product is a recombinant rabbit monoclonal antibody.


性能

形式	Liquid
存放说明	Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C.
解离常数 (K_D)	K _D = 5.80 x 10 ⁻¹¹ M



[Learn more about K_D](#)

存储溶液	pH: 7.20 Preservative: 0.05% Sodium azide Constituents: 0.1% BSA, 40% Glycerol, 9.85% Tris glycine, 50% Tissue culture supernatant
-------------	--



Western blot - Anti-Leukotriene A4 hydrolase antibody [EPR5712] (ab109434)

All lanes : Anti-Leukotriene A4 hydrolase antibody [EPR5712] (ab109434) at 1/10000 dilution

Lane 1 : HeLa cell lysate


Lane 2 : 293T cell lysate

Lane 3 : T47D cell lysate

Lane 4 : A549 cell lysate


Lysates/proteins at 10 µg per lane.

Predicted band size: 69 kDa



Immunohistochemical analysis of Leukotriene A4 hydrolase in paraffin-embedded Human breast carcinoma tissue using ab109434 at 1/250 dilution.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Leukotriene A4 hydrolase antibody [EPR5712] (ab109434)



Other - Anti-Leukotriene A4 hydrolase antibody [EPR5712] (ab109434)

Equilibrium dissociation constant (K_D)

Learn more about K_D

[Click here to learn more about \$K_D\$](#)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors