

Product datasheet

Anti-M6PR (cation independent) antibody [MEM-238] (Allophycocyanin) ab77023

概述

产品名称	Anti-M6PR (cation independent)抗体[MEM-238] (Allophycocyanin)
描述	小鼠单克隆抗体[MEM-238] to M6PR (cation independent) (Allophycocyanin)
宿主	Mouse
偶联物	Allophycocyanin. Ex: 645nm, Em: 660nm
经测试应用	适用于: Flow Cyt
种属反应性	与反应: Human 预测可用于: Non human primates
免疫原	Recombinant Vaccinia virus encoding CD222
阳性对照	Human blood cells.

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 15mM Sodium Azide Constituents: 0.2% BSA (high-grade protease free), PBS
纯度	Size exclusion
纯化说明	ab77023 is purified by size exclusion chromatography and adjusted for direct use.
克隆	单克隆
克隆编号	MEM-238
同种型	IgG1

应用

Our [Abpromise guarantee](#) covers the use of **ab77023** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
Flow Cyt		
应用说明	Flow Cyt: Use 10µl for 10 ⁶ cells or 100µl of whole blood. Not yet tested in other applications. Optimal dilutions/concentrations should be determined by the end user.	
靶标		
功能	Transport of phosphorylated lysosomal enzymes from the Golgi complex and the cell surface to lysosomes. Lysosomal enzymes bearing phosphomannosyl residues bind specifically to mannose-6-phosphate receptors in the Golgi apparatus and the resulting receptor-ligand complex is transported to an acidic prelysosomal compartment where the low pH mediates the dissociation of the complex. This receptor also binds IGF2. Acts as a positive regulator of T-cell coactivation, by binding DPP4.	
序列相似性	Belongs to the MRL1/IGF2R family. Contains 1 fibronectin type-II domain.	
结构域	Contains 15 repeating units of approximately 147 AA harboring four disulfide bonds each. The most highly conserved region within the repeat consists of a stretch of 13 AA that contains cysteines at both ends.	
细胞定位	Lysosome membrane. Colocalized with DPP4 in internalized cytoplasmic vesicles adjacent to the cell surface.	

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors