

Product datasheet

Anti-MAN2B1 antibody ab104521

2 图像

概述

产品名称	Anti-MAN2B1抗体
描述	兔多克隆抗体to MAN2B1
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Full length protein, amino acids 1 to 1011 of Human MAN2B1 (NP_000519.2).
阳性对照	MAN2B1 transfected 293T cell line lysate. HepG2 cell line lysate.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	pH: 7.20 Constituent: 2.68% PBS
纯度	Protein A purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab104521** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/500. Predicted molecular weight: 114 kDa.

靶标


相关性 MAN2B1 is necessary for the catabolism of N-linked carbohydrates released during glycoprotein

turnover. It cleaves all known types of alpha-mannosidic linkages.

细胞定位

Lysosome


图片



Western blot - Anti-MAN2B1 antibody (ab104521)

Anti-MAN2B1 antibody (ab104521) at 1/500 dilution + HepG2 cell line lysate at 50 µg

Predicted band size: 114 kDa



Western blot - Anti-MAN2B1 antibody (ab104521)

All lanes : Anti-MAN2B1 antibody (ab104521) at 1/500 dilution

Lane 1 : MAN2B1 transfected 293T cell line lysate

Lane 2 : Non transfected 293T cell line lysate

Lysates/proteins at 25 µg per lane.

Predicted band size: 114 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors