

Product datasheet

Anti-MAX binding protein antibody [MNTD1C3] ab53487

2 References

概述

产品名称	Anti-MAX binding蛋白抗体[MNTD1C3]
描述	小鼠单克隆抗体[MNTD1C3] to MAX binding蛋白
宿主	Mouse
经测试应用	适用于: WB, Dot blot
种属反应性	与反应: Human
免疫原	Recombinant fragment, internal sequence (Human).

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.05% Sodium Azide Constituents: 1% BSA, PBS, pH 7.4
纯度	Protein G purified
纯化说明	Purified by Protein G column from culture supernatant of hybridoma, cultured in medium containing bovine IgG depleted (approx. 95%) fetal bovine serum. Filtered through a 0.22µm membrane.
克隆	单克隆
克隆编号	MNTD1C3
同种型	IgG1

应用

Our [Abpromise guarantee](#) covers the use of **ab53487** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use at an assay dependent dilution. Detects a band of approximately 38 kDa (predicted molecular weight: 62 kDa).

应用	Ab评论	说明
----	------	----

Dot blot		Use at an assay dependent dilution.
----------	--	-------------------------------------

靶标

功能	Binds DNA as a heterodimer with MAX and represses transcription. Binds to the canonical E box sequence 5'-CACGTG-3' and, with higher affinity, to 5'-CACGCG-3'.
序列相似性	Contains 1 basic helix-loop-helix (bHLH) domain.
细胞定位	Nucleus.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors