

Product datasheet

Anti-MICAL1 antibody ab101529

1 图像

概述

产品名称	Anti-MICAL1 抗体
描述	兔多克隆抗体 to MICAL1
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Human 预测可用于: Mouse, Rat, Cow
免疫原	Recombinant fragment corresponding to a region within amino acids 225-401 of Human MICAL1 (AAH09972).
阳性对照	HCT116 whole cell lysate

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.01% Thimerosal (merthiolate) Constituents: 20% Glycerol, 0.1M Tris, 0.1M Glycine, pH 7.0
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab101529** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/500 - 1/3000. Predicted molecular weight: 109 kDa.

图标

功能	May be a cytoskeletal regulator that connects NEDD9 to intermediate filaments.
组织特异性	Expressed in the thymus, lung, spleen, kidney, testis and hematopoietic cells.
序列相似性	Contains 1 CH (calponin-homology) domain. Contains 1 LIM zinc-binding domain.
结构域	The C-terminal coiled coil part contains the plexin-interacting region.
细胞定位	Cytoplasm. Cytoplasm > cytoskeleton.

图片



Anti-MICAL1 antibody (ab101529) at 1/1000 dilution + HCT116 whole cell lysate at 30 µg

Predicted band size: 109 kDa

7.5% SDS-PAGE

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors