

Product datasheet

Anti-MLLT1 antibody ab71477

1 图像

概述

产品名称	Anti-MLLT1抗体
描述	兔多克隆抗体to MLLT1
宿主	Rabbit
经测试应用	适用于: ELISA, WB
种属反应性	与反应: Human
免疫原	Synthetic peptide from the C-terminal region of human MLLT1 conjugated to KLH
阳性对照	A257 cell line lysate.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: 0.09% Sodium Azide Constituents: PBS
纯度	Protein G purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab71477** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ELISA		1/1000.
WB		1/100 - 1/500. Detects a band of approximately 62 kDa (predicted molecular weight: 62 kDa).

靶标

功能	Capable of activating transcription from synthetic reporter genes in both lymphoid and myeloid cells.
疾病相关	Note=A chromosomal aberration involving MLLT1 is associated with acute leukemias. Translocation t(11;19)(q23;p13.3) with MLL/HRX. The result is a rogue activator protein.
序列相似性	Contains 1 YEATS domain.
细胞定位	Nucleus.

图片



Anti-MLLT1 antibody (ab71477) at 1/100 dilution + A257 cell line lysate at 35 µg

Predicted band size: 62 kDa

Observed band size: 62 kDa

Additional bands at: 75 kDa. We are unsure as to the identity of these extra bands.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors