

Product datasheet

Anti-MTHFD2L antibody ab90887

1 图像

概述

产品名称	Anti-MTHFD2L 抗体
描述	兔多克隆抗体 to MTHFD2L
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Human 预测可用于: Mouse, Rat, Rabbit, Horse, Chicken, Guinea pig, Cow, Cat, Dog, Pig, Zebrafish
免疫原	Synthetic peptide corresponding to a region within amino acids 216-265 (SKNVGMPIAM LLHTDGEHER PGGDATVTIA HRYTPKEQLK IHTQLADIII) of Human MTHFD2L (NP_001138450). Run BLAST with ExPASy Run BLAST with NCBI
阳性对照	Human fetal lung lysate

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: 0.09% Sodium azide Constituents: 2% Sucrose, PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab90887** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 31 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

靶标

组织特异性	Isoform 1, isoform 4 and isoform 5 are expressed in brain and placenta.
通路	One-carbon metabolism; tetrahydrofolate interconversion.
序列相似性	Belongs to the tetrahydrofolate dehydrogenase/cyclohydrolase family.
细胞定位	Mitochondrion inner membrane.

图片



Anti-MTHFD2L antibody (ab90887) at 1 µg/ml (in 5% skim milk / PBS buffer) + Human fetal lung lysate at 10 µg

Secondary

HRP conjugated anti-Rabbit IgG at 1/50000 dilution

Predicted band size: 31 kDa

Gel concentration: 12%

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors