

## Product datasheet

# Anti-Macrophage Specific Gene antibody ab25146

## 2 References

### 概述

产品名称	Anti-Macrophage Specific Gene抗体
描述	兔多克隆抗体to Macrophage Specific Gene
宿主	Rabbit
特异性	ab25146 recognises Mpg-1 (macrophage-specific gene, MSG, MPS1).
经测试应用	适用于: WB
种属反应性	与反应: Mouse, Rat
免疫原	Recombinant full length protein (Rat)
阳性对照	LPS-stimulated RAW264.7 and rat spleen.

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
存储溶液	Preservative: 0.1% Sodium Azide Constituents: PBS
纯度	Protein A purified
克隆	多克隆
同种型	IgG

### 应用

Our [Abpromise guarantee](#) covers the use of **ab25146** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/2000. Predicted molecular weight: 71-74 kDa.

### 图标

---

组织特异性	Expressed in macrophages and peripheral blood monocytes.
序列相似性	Contains 1 MACPF domain.
细胞定位	Membrane.

---

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

---

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

---

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors