

Product datasheet

Anti-NCKAP1L antibody ab82664

1 Abreviews 1 图像

概述

产品名称	Anti-NCKAP1L抗体
描述	兔多克隆抗体to NCKAP1L
宿主	Rabbit
经测试应用	适用于: WB, ELISA
种属反应性	与反应: Human 预测可用于: Mouse, Rat, Rabbit, Horse, Guinea pig, Cow, Cat, Dog
免疫原	A synthetic peptide corresponding to a region within the C terminal amino acids 1008-1057 (LPLLATDPSS FYSIEKDGYN NNIHCLTKAI IQVSAALFTL YNKNIETHLK) of human NCKAP1L (NP_005328). Run BLAST with ExPASy Run BLAST with NCBI
阳性对照	293T cell lysate

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
存储溶液	Preservative: 0.09% Sodium azide Constituents: 2% Sucrose, PBS
纯度	Immunogen affinity purified
纯化说明	Purified by peptide affinity chromatography method.
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab82664** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 128 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.
ELISA		Use at an assay dependent concentration. ELISA titre using peptide based assay, 1:62500.

靶标

组织特异性	Expressed only in cells of hematopoietic origin.
序列相似性	Belongs to the HEM-1/HEM-2 family.
细胞定位	Cell membrane.

图片



Anti-NCKAP1L antibody (ab82664) at 1 µg/ml
(in 5% skim milk / PBS buffer) + 293T cell
lysate at 10 µg

Secondary

HRP conjugated anti-Rabbit IgG diluted in
1/50000 - 1/100000 as secondary antibody.

Predicted band size: 128 kDa

Observed band size: 128 kDa

Gel concentration: 12%

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Promise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors