

## Product datasheet

# Anti-NDUFA4L2 antibody ab74138

★★★★☆ 2 Abreviews 3 References 1 图像

### 概述

产品名称	Anti-NDUFA4L2抗体
描述	兔多克隆抗体to NDUFA4L2
宿主	Rabbit
经测试应用	适用于: WB, ELISA
种属反应性	与反应: Mouse, Human
免疫原	Synthetic peptide derived from an internal sequence of human NDUFA4L2.
阳性对照	Extracts from Jurkat cells, COLO cells, HUVEC cells and HepG2 cells.

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	pH: 7.40 Preservative: 0.02% Sodium azide Constituents: 50% Glycerol, 0.87% Sodium chloride, PBS  Without Mg <sup>2+</sup> and Ca <sup>2+</sup>
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

### 应用

Our [Abpromise guarantee](#) covers the use of **ab74138** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明
WB	★★★★☆	1/500 - 1/1000. Detects a band of approximately <19 kDa (predicted molecular weight: 10 kDa).
ELISA		1/20000.

## 靶标

### 序列相似性

Belongs to the complex I NDUFA4 subunit family.

## 图片



**All lanes :** Anti-NDUFA4L2 antibody (ab74138) at 1/500 dilution

**Lane 1 :** Extracts from Jurkat cells

**Lane 2 :** Extracts from COLO cells

**Lane 3 :** Extracts from HUVEC cells

**Lane 4 :** Extracts from HepG2 cells

**Lane 5 :** Extracts from HepG2 cells with immunising peptide at 5 µg

Lysates/proteins at 5 µg per lane.

**Predicted band size:** 10 kDa

**Observed band size:** <19 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors