

Product datasheet

Anti-NEUROD4 antibody ab90484

1 图像

概述

产品名称	Anti-NEUROD4抗体
描述	兔多克隆抗体to NEUROD4
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Human 预测可用于: Mouse, Rat, Sheep, Horse, Chicken, Guinea pig, Cow, Cat, Dog, Pig, Zebrafish



**免疫原** Synthetic peptide designed within residues: HMETHLLHLK PQVFKSLGES SFGSHLPDCS TPPYEGPLTP PLSISGNFSL, corresponding to internal amino acids 216 - 265 of Human NEUROD4 (NP\_067014).

 [Run BLAST with ExPASy](#)  [Run BLAST with NCBI](#)

**阳性对照** Human fetal spleen lysate

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: 0.09% Sodium azide Constituents: PBS, 2% Sucrose
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab90484** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 37 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

## 靶标

相关性	Appears to mediate neuronal differentiation
细胞定位	Nuclear

## 图片



Anti-NEUROD4 antibody (ab90484) at 1 µg/ml (in 5% skim milk / PBS buffer) + Human fetal spleen lysate at 10 µg

### Secondary

HRP conjugated anti-Rabbit IgG at 1/50000 dilution

**Predicted band size:** 37 kDa

Gel concentration 12%

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors