

Product datasheet

Anti-NLRP10 antibody ab115966

1 图像

概述

产品名称	Anti-NLRP10抗体
描述	兔多克隆抗体to NLRP10
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	KLH conjugated synthetic peptide designed between amino acids 465-494 from the C-terminal region of Human NLRP10 (NP_789791.1).
阳性对照	Human placenta tissue lysates.

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at 4°C (up to 6 months). Store at -20°C long term.
存储溶液	Preservative: 0.09% Sodium azide Constituent: 99% PBS
纯度	Immunogen affinity purified
纯化说明	This antibody is purified through a protein A column, followed by peptide affinity purification.
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab115966** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/100 - 1/500. Predicted molecular weight: 75 kDa.

图标

<b>功能</b>	Inhibits CASP1/caspase-1-dependent IL1B secretion and PYCARD-mediated apoptosis, but not that induced by FAS or BID.
<b>组织特异性</b>	Widely expressed with highest levels in heart, brain and skeletal muscle.
<b>序列相似性</b>	Belongs to the NLRP family. Contains 1 DAPIN domain. Contains 1 NACHT domain.

## 图片



Anti-NLRP10 antibody (ab115966) at 1/100 dilution + Human placenta tissue lysates at 35 µg

**Predicted band size: 75 kDa**

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors