

# Anti-NK-3R antibody ab123303

[1 References](#) [1 图像](#)

### 概述

|       |  |
|-------|--|
| 产品名称  | Anti-NK-3R抗体   |
| 描述    | 兔多克隆抗体to NK-3R   |
| 宿主    | Rabbit   |
| 经测试应用 | 适用于: WB  |
| 种属反应性 | 与反应: Human   |
| 免疫原   | Synthetic peptide corresponding to Human NK-3R aa 405-435 (C terminal) conjugated to keyhole limpet haemocyanin.<br>(Peptide available as <a href="#">ab176476</a> ) |
| 阳性对照  | K562 cell line lysates.  |
| 常规说明  | Protein previously labeled as Neurokinin 3 Receptor NK3.   |

### 性能

|      |   |
|------|---|
| 形式   | Liquid  |
| 存放说明 | Shipped at 4°C. Store at 4°C (up to 6 months). Store at -20°C long term.                    |
| 存储溶液 | Preservative: 0.09% Sodium azide<br>Constituent: 99% PBS                                    |
| 纯度   | Immunogen affinity purified   |
| 纯化说明 | ab123303 is purified through a protein A column, followed by peptide affinity purification. |
| 克隆   | 多克隆   |
| 同种型  | IgG   |

### 应用

#### The Abpromise guarantee

**Abpromise™** 承诺保证使用ab123303于以下的经测试应用

“应用说明”部分下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

| 应用 | Ab评论 | 说明   |
|----|------|--|
| WB |      | 1/100 - 1/500. Predicted molecular weight: 52 kDa. |

## 靶标


---

|       |   |
|-------|---|
| 功能    | This is a receptor for the tachykinin neuropeptide neuromedin-K (neurokinin B). It is associated with G proteins that activate a phosphatidylinositol-calcium second messenger system. The rank order of affinity of this receptor to tachykinins is: neuromedin-K > substance K > substance P. |
| 疾病相关  | Hypogonadotropic hypogonadism 11 with or without anosmia  |
| 序列相似性 | Belongs to the G-protein coupled receptor 1 family.   |
| 翻译后修饰 | The anchoring of this receptor to the plasma membrane is probably mediated by the palmitoylation of a cysteine residue.   |
| 细胞定位  | Cell membrane.  |

---

## 图片

---



Anti-NK-3R antibody (ab123303) at 1/100 dilution + K562 cell line lysates at 35  $\mu$ g

**Predicted band size:** 52 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors