

## Product datasheet

# Anti-ODF3 antibody ab111496

### 1 图像

#### 概述

产品名称	Anti-ODF3抗体
描述	兔多克隆抗体to ODF3
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Mouse 预测可用于: Human
免疫原	Synthetic peptide conjugated to KLH, corresponding to a region within C terminal amino acids 177-207 of Human ODF3 (NP_444510.2).
阳性对照	Mouse testis tissue lysate

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
存储溶液	Preservative: 0.09% Sodium azide Constituent: PBS
纯度	Immunogen affinity purified
纯化说明	ab111496 is purified through a protein A column, followed by peptide affinity purification.
克隆	多克隆
同种型	IgG

#### 应用

Our [Abpromise guarantee](#) covers the use of **ab111496** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/100 - 1/500. Predicted molecular weight: 28 kDa.

## 靶标

功能	Outer dense fibers are filamentous structures located on the outside of the axoneme in the midpiece and principal piece of the mammalian sperm tail. May help to maintain the passive elastic structures and elastic recoil of the sperm tail.
组织特异性	Testis-specific.
序列相似性	Belongs to the ODF3 family. Contains 2 DUF1309 repeats.
细胞定位	Cytoplasm. Expressed in the cytoplasmic lobe of spermatids.

## 图片



Anti-ODF3 antibody (ab111496) at 1/100 dilution + Mouse testis tissue lysate at 35 µg

**Predicted band size:** 28 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors