

Product datasheet

Anti-PCDHB10 antibody ab121400

1 图像

概述

产品名称	Anti-PCDHB10抗体
描述	兔多克隆抗体to PCDHB10
宿主	Rabbit
经测试应用	适用于: IHC-P
种属反应性	与反应: Human
免疫原	antigen, corresponding to amino acids 45-177 of Human PCDHB10.
阳性对照	Human nasopharynx tissue.

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
存储溶液	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 59% PBS, 40% Glycerol
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab121400** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab 评论	说明
IHC-P		1/50 - 1/200. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol. Perform Heat Induced Epitope Retrieval (HIER) at pH 6 using a pressure boiler as heat source.

## 靶标


---

功能	Potential calcium-dependent cell-adhesion protein. May be involved in the establishment and maintenance of specific neuronal connections in the brain.
序列相似性	Contains 6 cadherin domains.
细胞定位	Cell membrane.

---

## 图片

---



ab121400 at 1/150 dilution staining  
PCDHB10 in Paraffin-embedded Human  
nasopharynx tissue by Immunohistochemistry.  
Note: strong cytoplasmic positivity in  
respiratory epithelial cells.

Immunohistochemistry (Formalin/PFA-fixed paraffin-  
embedded sections) - Anti-PCDHB10 antibody  
(ab121400)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors