

Product datasheet

Anti-PCNP antibody ab97909

2 图像

概述

产品名称	Anti-PCNP抗体
描述	兔多克隆抗体to PCNP
宿主	Rabbit
经测试应用	适用于: WB, IHC-P
种属反应性	与反应: Human
免疫原	Synthetic peptide, corresponding to a sequence within amino acids 1-20 of Human PCNP (NP_065090). Run BLAST with ExPASy Run BLAST with NCBI
阳性对照	293T, H1299, HeLa, HepG2, MOLT4 and Raji whole cell lysates; H1299 xenograft tissue.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	pH: 7.00 Preservative: 0.01% Thimerosal (merthiolate) Constituents: 2.68% PBS, 10% Glycerol
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab97909** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/500 - 1/3000. Predicted molecular weight: 19 kDa.
IHC-P		1/100 - 1/500.

靶标

功能	May be involved in cell cycle regulation.
翻译后修饰	Ubiquitinated; mediated by UHRF2 and leading to its subsequent proteasomal degradation.
细胞定位	Nucleus.

图片



All lanes : Anti-PCNP antibody (ab97909) at 1/1000 dilution


Lane 1 : H1299 whole cell lysate

Lane 2 : HeLa whole cell lysate

Lysates/proteins at 30 µg per lane.

Predicted band size: 19 kDa

12% SDS-PAGE.



ab97909 at 1/500 dilution staining PCNP in H1299 xenograft by Immunohistochemistry, Paraffin-embedded tissue.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors