

Product datasheet

Anti-PDIA5 antibody ab97328

1 图像

概述

产品名称	Anti-PDIA5抗体
描述	兔多克隆抗体to PDIA5
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Human 预测可用于: Mouse, Rat, Cow 
免疫原	Recombinant fragment within Human PDIA5 aa 1-223. The exact sequence is proprietary. (AAH01625) Database link: <a href="#">Q14554</a>
阳性对照	WB: HepG2, A549 and HCT116 cell lysates.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	pH: 7.00 Preservative: 0.01% Thimerosal (merthiolate) Constituents: 0.75% Glycine, 20% Glycerol, 1.21% Tris
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab97328** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/500 - 1/3000. Predicted molecular weight: 60 kDa.

## 靶标

### 序列相似性


Belongs to the protein disulfide isomerase family.

Contains 3 thioredoxin domains.

### 细胞定位

Endoplasmic reticulum lumen.

## 图片



Anti-PDIA5 antibody (ab97328) at 1/1000 dilution + HepG2 whole cell lysate at 30 µg

**Predicted band size:** 60 kDa

12% SDS-PAGE.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors