

Product datasheet

PSME1 peptide ab4949

概述

产品名称 PSME1多肽

描述

性质 Synthetic

技术指标

Our [Abpromise guarantee](#) covers the use of **ab4949** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用 Blocking

形式 Liquid

补充说明

This peptide may be used for neutralization and control experiments with the polyclonal antibody that reacts with this product and human proteasome activator 11S REG alpha, catalog [ab3333](#). Using a solution of peptide of equal volume and concentration to the corresponding antibody will yield a large molar excess of peptide (~ 70-fold) for competitive inhibition of antibody-protein binding reactions.

制备和贮存

稳定性和存储 Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

常规信息

功能 Implicated in immunoproteasome assembly and required for efficient antigen processing. The PA28 activator complex enhances the generation of class I binding peptides by altering the cleavage pattern of the proteasome.

序列相似性 Belongs to the PA28 family.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors