

# PTIP overexpression 293T lysate (whole cell) ab94247

## 2 图像

### 概述

<b>产品名称</b>	PTIP overexpression 293T裂解物(whole cell)
<b>常规说明</b>	ab94247 is a 293T cell transfected lysate in which Human PTIP has been transiently over-expressed using a pCMV-PTIP plasmid. The lysate is provided in 1X Sample Buffer. <b>Note:</b> For more details about how the transfected lysate was prepared <a href="#">view preparation notes</a>
<b>经测试应用</b>	<b>适用于:</b> WB

### 性能

<b>形式</b>	Liquid
<b>存放说明</b>	Shipped on dry ice. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
<b>存储溶液</b>	pH: 6.80 Constituent: 100% 1x Sample Buffer
<b>纯度</b>	Whole Cell Lysate
<b>背景</b>	<p>Function: Involved in DNA damage response and in transcriptional regulation through histone methyltransferase (HMT) complexes. Plays a role in early development. In DNA damage response is required for cell survival after ionizing radiation. In vitro shown to be involved in the homologous recombination mechanism for the repair of double-strand breaks (DSBs). Its localization to DNA damage foci requires RNF8 and UBE2N. Recruits TP53BP1 to DNA damage foci and, at least in particular repair processes, effective DNA damage response appears to require the association with TP53BP1 phosphorylated by ATM at 'Ser-25'. Together with TP53BP1 regulates ATM association. Recruits PA1 to sites of DNA damage and the PA1:PAXIP1 complex is required for cell survival in response to DNA damage; the function is probably independent of MLL-containing histone methyltransferase (HMT) complexes. Promotes ubiquitination of PCNA following UV irradiation and may regulate recruitment of polymerase eta and RAD51 to chromatin after DNA damage. Proposed to be involved in transcriptional regulation by linking MLL-containing histone methyltransferase (HMT) complexes to gene promoters by interacting with promoter-bound transcription factors such as PAX2. Associates with gene promoters that are known to be regulated by MLL2. During immunoglobulin class switching in activated B cells is involved in trimethylation of histone H3 at 'Lys-4' and in transcription initiation of downstream switch regions at the immunoglobulin heavy-chain (Igh) locus; this function appears to involve the recruitment of MLL-containing HMT complexes.</p> <p>Similarity: Contains 6 BRCT domains. Domain: The BRCT 5 and 6 domains function as a single module and are necessary and sufficient for in vitro phospho-specific binding (substrates phosphorylated by the kinases ataxia telangiectasia-mutated (ATM), ataxia telangiectasia and RAD3-related (ATR) in response to gamma irradiation). In contrast, in vivo two pairs of BRCT</p>

domains (3-6) bind to phosphorylated TP53BP1 much more efficiently.


## 应用

Our [Abpromise guarantee](#) covers the use of **ab94247** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use at an assay dependent dilution.

## 图片



ab94247 at 15µg/lane on an SDS-PAGE gel.



**All lanes** : Anti-PTIP antibody ([ab56934](#)) at  
1/500 dilution

**Lane 1** : PTIP overexpression 293T lysate  
(whole cell) (ab94247)

**Lane 2** : 293T non-transfected lysate

Lysates/proteins at 25 µg per lane.

### Secondary

**All lanes** : Goat Anti-mouse IgG (H and L)  
HRP conjugated at 1/2500 dilution

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet

- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors