

Product datasheet

Anti-Peripherin antibody ab123576

1 图像

概述

产品名称	Anti-Peripherin抗体
描述	兔多克隆抗体to Peripherin
宿主	Rabbit
经测试应用	适用于: WB, ICC/IF
种属反应性	与反应: Mouse, Rat, Cat, Human
免疫原	Recombinant full length Peripherin, purified from E. coli.
阳性对照	Mixed neuron/glia cultures from newborn Rat brain; whole Rat brain stem homogenate.

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Constituent: 99% Whole serum
纯度	Whole antiserum
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab123576** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab 评论	说明
WB		1/20000. Predicted molecular weight: 54 kDa.
ICC/IF		1/1000 - 1/5000. 1:1000 using fluorescent secondary antibodies or 1:5000 using peroxidase or other enzyme linked methods immunocytochemical procedures.

## 靶标

---


**功能** Class-III neuronal intermediate filament protein.

**序列相似性** Belongs to the intermediate filament family.

---

## 图片

---



ab123576, at 1/1000, staining Peripherin (green) in a mixed neuron/glia culture from newborn Rat brain by immunofluorescence. Red shows staining for phosphorylated NF-H, while blue indicates the localization of DNA.

Immunocytochemistry/ Immunofluorescence - Anti-Peripherin antibody (ab123576)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors