

Anti-RNF186 antibody ab86547

[1 References](#) [1 图像](#)

概述

产品名称	Anti-RNF186抗体
描述	兔多克隆抗体to RNF186
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Human 预测可用于: Horse, Guinea pig 
免疫原	Synthetic peptide corresponding to a region within the internal amino acids 108-157 (QLAQPCTEVS LCPQGLVDPA DLAAGHPSLV GEDGQDEVSA NHVAARRLAA) of human RNF186 (NP_061935). Run BLAST with ExPASy Run BLAST with NCBI
阳性对照	WB: Transfected 293T cell lysate.

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
存储溶液	pH: 7.2 Preservative: 0.09% Sodium azide Constituents: PBS, 2% Sucrose
纯度	Immunogen affinity purified
纯化说明	Purified by peptide affinity chromatography method.
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab86547于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。


应用	Ab评论	说明
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 24 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

靶标

相关性 RNF186 contains a RING type zinc finger, a motif known to be involved in protein-protein and protein-DNA interactions. The specific function of this protein is unknown.

细胞定位 Membrane; Multi pass membrane protein.

图片



Western blot - Anti-RNF186 antibody (ab86547)

Anti-RNF186 antibody (ab86547) at 1 µg/ml (in 5% skim milk / PBS buffer) + transfected 293T cell lysate at 10 µg

Secondary

HRP conjugated anti-Rabbit IgG at 1/50000 dilution

Predicted band size: 24 kDa

Observed band size: 24 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors