

Product datasheet

Rabbit Anti-Chicken IgY H&L ab31170

1 References 1 图像

概述

产品名称	兔抗鸡IgY H&L
宿主	Rabbit
靶标种属	Chicken
特异性	ab31170 recognises chicken IgY.
经测试应用	<b>适用于:</b> Immunoelectrophoresis, WB, Immunodiffusion, ELISA
免疫原	Purified chicken egg yolk IgY (H/L)

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 0.0268% PBS, 0.812% Sodium chloride
纯度	Immunogen affinity purified
纯化说明	Determined by SDS-PAGE(blue stain), showing that the IgG constitute was >90% of total protein.
克隆	多克隆
同种型	IgG

应用


Our [Abpromise guarantee](#) covers the use of **ab31170** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
Immunoelectrophoresis		Use at an assay dependent dilution.
WB		Use at an assay dependent dilution. Predicted molecular weight: 10 kDa.
Immunodiffusion		Use at an assay dependent dilution.

应用	Ab评论	说明
ELISA		1/200000.

## 图片



Western blot - Rabbit polyclonal Secondary  
Antibody to Chicken IgY - H&L (ab31170)

Ab31170 was used at dilution 1/4000 with the primary antibody [ab26184](#) in WB. See the review on [ab26184](#).

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors