

## Product datasheet

# Rabbit Anti-Human IgG H&L (Alkaline Phosphatase) ab97157

### 概述

产品名称	兔抗人IgG H&L (Alkaline Phosphatase)
宿主	Rabbit
靶标种属	Human
特异性	By immunoelectrophoresis and ELISA this antibody reacts specifically with Human IgG and with light chains common to other Human immunoglobulins. No antibody was detected against non immunoglobulin serum proteins.
经测试应用	适用于: ICC, IHC-P, ELISA, WB
偶联物	Alkaline Phosphatase

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 7.10 Preservative: 0.09% Sodium azide Constituents: 0.2% BSA, 0.58% Sodium chloride, 0.0095% Magnesium chloride, 0.00136% Zinc chloride, 1.19% HEPES
纯度	Immunogen affinity purified
纯化说明	This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to Alkaline Phosphatase (AP).
克隆	多克隆
同种型	IgG

### 应用

Our [Abpromise guarantee](#) covers the use of **ab97157** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ICC		Use at an assay dependent dilution.
IHC-P		1/100 - 1/1000.

应用	Ab评论	说明
ELISA		1/5000 - 1/50000. (as Primary).
WB		1/5000 - 1/50000. (Chemiluminescent).

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors