

Product datasheet

Recombinant Human DYNLRB1 protein ab98120

1 图像

概述

产品名称	重组人DYNLRB1蛋白
蛋白长度	Full length protein

描述

性质	Recombinant
来源	Escherichia coli

氨基酸序列

Accession	<a href="#">Q9NP97</a>
种属	Human
序列	MAEVEETLKR LQSQKGVQGI IVVNTEGIPI KSTMDNPTTT QYASLMHSFI LKARSTVRDI DPQNDLTFLR IRSKKNEIMV APDKDYFLIVIQNPTLEHH HHHH

分子量	12 kDa including tags
氨基酸	1 to 96
标签	His tag C-Terminus

技术指标

Our [Abpromise guarantee](#) covers the use of **ab98120** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE Mass Spectrometry
质谱法	MALDI-TOF
纯度	> 95 % SDS-PAGE. ab98120 is purified by using conventional chromatography techniques.
形式	Liquid

制备和贮存

## 稳定性和存储

Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

pH: 8.00

Constituents: 0.316% Tris HCl, 10% Glycerol

## 常规信息

### 相关性

DYNLRB1 is a member of the roadblock dynein light chain family and encodes a cytoplasmic protein that is capable of binding intermediate chain proteins. Upregulation of this gene has been associated with hepatocellular carcinomas, suggesting that this gene may be involved in tumor progression.

### 细胞定位

Cytoplasmic

## 图片



15% SDS-PAGE analysis of ab98120 (3ug)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors