

Product datasheet

Recombinant Human Dynein light chain protein ab107140

1 图像

概述

产品名称	重组人Dynein light chain蛋白
蛋白长度	Full length protein

描述

性质	Recombinant
来源	Escherichia coli
氨基酸序列	
Accession	O96015
种属	Human
序列	<p>MGSSHHHHHSSGLVPRGSHMGETEGKKDEADYKRLQT</p> <p>FPLVRHSDMPEE</p> <p>MRVETMELCVTACEKFSNNNESAAKMIKETMDKKFGSSWH</p> <p>VVIGEGFGFE ITHEVKNLLYLYFGGTLAVCVWKCS</p>
分子量	14 kDa including tags
氨基酸	1 to 105
标签	His tag N-Terminus

技术指标

Our [Abpromise guarantee](#) covers the use of **ab107140** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	<p>SDS-PAGE</p> <p>Mass Spectrometry</p>
质谱法	MALDI-TOF
纯度	<p>> 95 % SDS-PAGE.</p> <p>ab107140 is purified using conventional chromatography techniques.</p>
形式	Liquid

制备和贮存

稳定性和存储

Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.


pH: 8.00

Constituents: 0.0154% DTT, 0.316% Tris HCl, 20% Glycerol, 0.58% Sodium chloride

常规信息

功能	Force generating protein of respiratory cilia. Produces force towards the minus ends of microtubules. Dynein has ATPase activity.
序列相似性	Belongs to the dynein light chain family.
细胞定位	Cytoplasm > cytoskeleton > cilium axoneme.

图片



15% SDS-PAGE showing ab107140 (3 µg) at approximately 14.1 kDa.

SDS-PAGE - Recombinant Human Dynein light chain protein (ab107140)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors