

## Product datasheet

# Recombinant Human EBF1 protein ab112729

### 概述

---

产品名称	重组人EBF1蛋白
蛋白长度	Protein fragment

### 描述

---

性质	Recombinant
来源	Escherichia coli

### 氨基酸序列

Accession	<a href="#">Q9UH73</a>
种属	Human
分子量	36 kDa

### 技术指标

---

Our [Abpromise guarantee](#) covers the use of **ab112729** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Western blot
形式	Liquid

### 制备和贮存

---

稳定性和存储	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles. pH: 6.80 Constituents: 1.58% Tris HCl, 3.09% DTT, 4% SDS, 0.2% Bromophenol blue, 20% Sucrose
--------	---

### 常规信息

---

功能	Transcriptional activator which recognizes variations of the palindromic sequence 5'-ATTCCCNNGGGAATT-3'.
序列相似性	Belongs to the COE family. Contains 1 IPT/TIG domain.

---

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

---

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

---

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors