

Product datasheet

Recombinant Human ELK4 protein ab114409

1 图像

概述

产品名称	重组人ELK4蛋白
蛋白长度	Protein fragment

描述

性质	Recombinant
来源	Wheat germ

氨基酸序列

Accession	P28324
种属	Human
序列	KDVENGGKDKPPQPGAKTSSRNDYIHSGLYSSFTLNSLNSSNVKLFKLIK TENPAEKLAEEKSPQEPTPSVIK FVTTPSKKPPVEPVAA
分子量	35 kDa
氨基酸	118 to 206

技术指标

Our [Abpromise guarantee](#) covers the use of **ab114409** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE ELISA Western blot
形式	Liquid
补充说明	Protein concentration is above or equal to 0.05 mg/ml.


制备和贮存

稳定性和存储	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.3% Glutathione, 0.79% Tris HCl
--------	---

常规信息

功能	Forms a ternary complex with the serum response factor (SRF). Requires DNA-bound SRF for ternary complex formation and makes extensive DNA contacts to the 5' side of SRF, but does not bind DNA autonomously.
序列相似性	Belongs to the ETS family. Contains 1 ETS DNA-binding domain.
细胞定位	Nucleus.

图片



12.5% SDS-PAGE analysis of ab114409 stained with Coomassie Blue.

SDS-PAGE - ELK4 protein (ab114409)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors