

Product datasheet

Recombinant Human GATA2 protein ab114360

1 图像

概述

产品名称	重组人GATA2蛋白
蛋白长度	Full length protein

描述

性质	Recombinant
来源	Wheat germ

氨基酸序列

Accession [P23769](#)

种属 Human

序列 MEVAPEQPRWMAHPAVLNAQHPDSHHPGLAHNYMEPAQLLPPDEVDFFN
 HLDSQGNPYYPANPAHARARVSYSPAARLTGGQMCRPHLLHSPGLPWLDG
 GKAALSAAAAHHHPWTVSPFSKTPLHPSAAGGPGGPLSVYPGAGGGSGG
 GSGSSVASLTPTAAHSGSHLFGFPPTPPKEVSPDPSTTGAASPASSSAGG
 SAARGEDKDGVKYQVSLTESMKMESGSPLRPLATMGTPATHHPIPTYP
 SYVPAAAHDYSSGLFHPGGFLGGPASSFTPQQRSKARSCSEGRECVNCGA
 TATPLWRRDGTGHYLCNACGLYHKMNGQNRPLIKPKRRLSAARRAGTCCA
 NCQTTTTTLWRRNANGDPVCNACGLYYKLHNVRPLTMKKEGIQTRNRKM
 SNKSKSKKGAECFEELSKCMQEKSSPFSAAALAGHMAPVGHLPFFSHSG
 HILPTPTPIHPSSSLSGHPPSSMVTAMG

分子量 79 kDa including tags

氨基酸 1 to 480

技术指标

Our [Abpromise guarantee](#) covers the use of **ab114360** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	ELISA SDS-PAGE Western blot
形式	Liquid

补充说明

Protein concentration is above or equal to 0.05 mg/ml. Best use within three months from the date of receipt of this protein.

制备和贮存

稳定性和存储

Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.

pH: 8.00

Constituents: 0.3% Glutathione, 0.79% Tris HCl

常规信息

功能

Transcriptional activator which regulates endothelin-1 gene expression in endothelial cells. Binds to the consensus sequence 5'-AGATAG-3'.

组织特异性

Endothelial cells.


序列相似性

Contains 2 GATA-type zinc fingers.

细胞定位

Nucleus.

图片



12.5% SDS-PAGE showing ab114360 at approximately 78.87kDa stained with Coomassie Blue.

SDS-PAGE - Recombinant Human GATA2 protein
(ab114360)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise,

please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors