

Product datasheet

Recombinant Human IFIT2 protein ab126008

概述

产品名称	重组人IFIT2蛋白
蛋白长度	Protein fragment

描述

性质	Recombinant
来源	Escherichia coli
氨基酸序列	
Accession	P09913
种属	Human
分子量	26 kDa
氨基酸	100 to 321

技术指标

Our [Abpromise guarantee](#) covers the use of **ab126008** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE
形式	Lyophilised

制备和贮存

稳定性和存储	Shipped at 4°C. Store at -20°C. Constituents: 0.32% Tris HCl, 0.58% Sodium chloride
复溶	Reconstitute with water to desired concentration.

常规信息

序列相似性	Belongs to the IFIT family. Contains 6 TPR repeats.
-------	--

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors