

Product datasheet

Recombinant Human Inhibin beta A protein ab53506

1 图像

概述

产品名称	重组人Inhibin beta A蛋白
蛋白长度	Full length protein

描述

性质	Recombinant
来源	Escherichia coli
氨基酸序列	
Accession	<a href="#">P08476</a>
种属	Human
氨基酸	311 to 426
标签	His tag N-Terminus
额外的序列信息	The sequence is the mature full length without the signal peptide and propeptide. Hexahistidine tag at N-terminus.

技术指标

Our [Abpromise guarantee](#) covers the use of **ab53506** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE
纯度	> 90 % SDS-PAGE.
形式	Liquid

制备和贮存

稳定性和存储	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
--------	---

常规信息

功能	Inhibins and activins inhibit and activate, respectively, the secretion of follitropin by the pituitary
----	---

gland. Inhibins/activins are involved in regulating a number of diverse functions such as hypothalamic and pituitary hormone secretion, gonadal hormone secretion, germ cell development and maturation, erythroid differentiation, insulin secretion, nerve cell survival, embryonic axial development or bone growth, depending on their subunit composition. Inhibins appear to oppose the functions of activins.


**序列相似性**

Belongs to the TGF-beta family.

**细胞定位**

Secreted.

**图片**



**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

**Our Abpromise to you: Quality guaranteed and expert technical support**

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

**Terms and conditions**

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors