

Product datasheet

Recombinant Human MAF1 protein ab111643

1 图像

概述

产品名称	重组人MAF1蛋白
蛋白长度	Full length protein

描述

性质	Recombinant
来源	Escherichia coli
氨基酸序列	
Accession	Q9H063
种属	Human

序列	<p>MGSSHHHHHHSSGLVPRGSHMKLLENSSF EAINSQLTV ETGDAHIIGRIE SYSCKMAGDDKHMFKQFCQEGQPHVLEALSPPQTSGLSPS RLSKSQGGEE EGPLSDKCSRKTLFYLIATLNESFRPDYDFSTARSH EFSR EPSLSWVVNA VNCSLFSAVREDFKDLKPQLWNAVDEEICLAECDIYSYNP DLDSDPFGED GSLWSFNFFYFNKRLKRIVFFSCRSISGSTYTPSEAGNEL DMELGEEVE EESRSGSGAEETSTMEEDRVPVICI</p>
----	--

分子量	31 kDa including tags
氨基酸	1 to 256
标签	His tag N-Terminus

技术指标

Our [Abpromise guarantee](#) covers the use of **ab111643** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Mass Spectrometry
	SDS-PAGE
质谱法	MALDI-TOF
纯度	> 90 % SDS-PAGE.

ab111643 is purified using conventional chromatography.

形式 Liquid

制备和贮存

稳定性和存储 Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 8.00

Constituent: 10% Glycerol

Note: Contains 90% Tris-HCl buffer.

常规信息


功能 Element of the mTORC1 signaling pathway that acts as a mediator of diverse signals and that represses RNA polymerase III transcription. Inhibits the de novo assembly of TFIIIB onto DNA.

序列相似性 Belongs to the MAF1 family.

翻译后修饰 Phosphorylated at Ser-60, Ser-68 and Ser-75; the major sites of phosphorylation. Nuclear accumulation correlates with a concomitant dephosphorylation. Phosphorylation may attenuate its RNA polymerase III-repressive function.

细胞定位 Nucleus. Cytoplasm.

图片



15% SDS-PAGE analysis of ab111643 (3 µg).

SDS-PAGE - Recombinant Human MAF1 protein
(ab111643)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors