

Product datasheet

Recombinant Human PAX9 protein ab152592

1 图像

概述

产品名称	重组人PAX9蛋白
蛋白长度	Full length protein

描述

性质	Recombinant
来源	Wheat germ

氨基酸序列

种属	Human
序列	MEPAFGEVNQLGGVFNVRPLPNAIRLRIVELAQLGIRPCDISRQLRVSH GCVSKILARYNETGSILPGAIGGSKPRVTTPTVVKHIRTYKQRDPGIFAW EIRDRLADGVCDKYNVPSVSSISRILRNKIGNLAQQGHYDSYKQHQP TPQALPYNHIYSYSPITAAAAKVPTPPGVPAIPGSVAMPRTWPSSHSVTD ILGIRSITDQVSDSSPYHSPKVEEWSSLGRNNFPAAAPHAVNGLEKGALE QEAKYQGAPNGLPAVGSFVSASSMAPYPTPAQVSPYMTYSAAPSGYVAGH GWQHAGGTSLSPHNCIPASLAFKGMQAAREGSHSVTASAL

氨基酸	1 to 341
标签	proprietary tag N-Terminus

技术指标

Our [Abpromise guarantee](#) covers the use of **ab152592** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Western blot ELISA
形式	Liquid
补充说明	Protein concentration is above or equal to 0.05 mg/ml.

制备和贮存

稳定性和存储	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.
--------	---


pH: 8.00

Constituents: 0.31% Glutathione, 0.79% Tris HCl

常规信息

功能	Transcription factor required for normal development of thymus, parathyroid glands, ultimobranchial bodies, teeth, skeletal elements of skull and larynx as well as distal limbs.
疾病相关	Defects in PAX9 are the cause of tooth agenesis selective type 3 (STHAG3) [MIM:604625]. A form of selective tooth agenesis, a common anomaly characterized by the congenital absence of one or more teeth. Selective tooth agenesis without associated systemic disorders has sometimes been divided into 2 types: oligodontia, defined as agenesis of 6 or more permanent teeth, and hypodontia, defined as agenesis of less than 6 teeth. The number in both cases does not include absence of third molars (wisdom teeth).
序列相似性	Contains 1 paired domain.
细胞定位	Nucleus.

图片



ab152592 on a 12.5% SDS-PAGE stained with Coomassie Blue.

SDS-PAGE - Recombinant Human PAX9 protein (ab152592)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors