

Recombinant Human Syndecan 4 protein ab117191

描述

产品名称	重组人Syndecan 4蛋白
纯度	> 85 % SDS-PAGE. ab117191 is purified by proprietary chromatographic techniques and is greater than 85% as determined by SDS-PAGE.
表达系统	Escherichia coli
Accession	<u>P31431</u>
蛋白长度	Protein fragment
无动物成分	No
性质	Recombinant
种属	Human
序列	MESIRETEVIDPQDLLEGYFSGALPDDEDVVGGQESDDFE LSGSGDLD DLED SMIGPEVVHPLVPLDNHIPERAGSGSQVPTEPKKLEEN EVIPKRIS PVEESEDVSNKVSMSSSTVQGSNIFERTE
预测分子量	14 kDa
氨基酸	19 to 145

技术指标

Our **Abpromise guarantee** covers the use of **ab117191** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE
形式	Liquid

制备和贮存

稳定性和存储	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle. pH: 8.00 Constituents: 10% Glycerol, 0.24% Tris
--------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

常规信息

功能	Cell surface proteoglycan that bears heparan sulfate.
组织特异性	Expressed in epithelial and fibroblastic cells.
序列相似性	Belongs to the syndecan proteoglycan family.
翻译后修饰	Shedding is enhanced by a number of factors such as heparanase, thrombin or EGF. Also by stress and wound healing. PMA-mediated shedding is inhibited by TIMP3.
细胞定位	Membrane. Secreted. Shedding of the ectodomain produces a soluble form.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors