

Product datasheet

Recombinant Human Tbx3 protein ab114501

1 图像

概述

产品名称	重组人Tbx3蛋白
蛋白长度	Protein fragment

描述

性质	Recombinant
来源	Wheat germ
氨基酸序列	
Accession	O15119
种属	Human
序列	KENGTSDESSSEQAAFNCFAQASSPAASTVGTSLKDLCP SEGESDAEAE SKEEHGPEACDAAKISTTTSEEPCKDKGSPAVKAHLFAAE RPRDSGRLDK
分子量	37 kDa including tags
氨基酸	331 to 430

技术指标

Our [Abpromise guarantee](#) covers the use of **ab114501** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	ELISA SDS-PAGE Western blot
形式	Liquid
补充说明	Protein concentration is above or equal to 0.05 mg/ml. This protein is best used within three months from the date of receipt.

制备和贮存

稳定性和存储	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.
--------	---


pH: 8.00

Constituents: 0.3% Glycerol, 0.79% Tris HCl

常规信息

功能	Transcriptional repressor involved in developmental processes. Probably plays a role in limb pattern formation.
组织特异性	Widely expressed.
疾病相关	Defects in TBX3 are the cause of ulnar-mammary syndrome (UMS) [MIM:181450]. UMS is characterized by ulnar ray defects, obesity, hypogonadism, delayed puberty, hypoplasia of nipples and apocrine glands.
序列相似性	Contains 1 T-box DNA-binding domain.
细胞定位	Nucleus.

图片



ab114501 analysed by 12.5% SDS-PAGE and stained with Coomassie Blue.

SDS-PAGE - Recombinant Human Tbx3 protein
(ab114501)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors