

Product datasheet

Recombinant Rat Choline Acetyltransferase protein ab119467

概述

| | |
|------|---------------------------------|
| 产品名称 | 重组大鼠Choline Acetyltransferase蛋白 |
| 蛋白长度 | Full length protein |

描述

| | |
|----|------------------|
| 性质 | Recombinant |
| 来源 | Escherichia coli |

氨基酸序列

| | |
|-----------|------------------------|
| Accession | P32738 |
| 种属 | Rat |
| 分子量 | 72 kDa |
| 氨基酸 | 1 to 640 |

技术指标

Our [Abpromise guarantee](#) covers the use of **ab119467** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| | |
|----|--|
| 应用 | Immunohistochemistry (Frozen sections) Blocking |
| 形式 | Liquid |

制备和贮存

| | |
|--------|---|
| 稳定性和存储 | Shipped at 4°C. Store at 4°C (up to 6 months). Store at -20°C. Constituents: 89% PBS, 10% Glycerol |
|--------|---|

常规信息

| | |
|------|--|
| 功能 | Catalyzes the reversible synthesis of acetylcholine (ACh) from acetyl CoA and choline at cholinergic synapses. |
| 疾病相关 | Myasthenic syndrome, congenital, 6, presynaptic |

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors