

Product datasheet

Recombinant human RAB37 protein ab90792

概述

产品名称	重组人RAB37蛋白
蛋白长度	Full length protein

描述

性质	Recombinant
来源	Escherichia coli

氨基酸序列

种属	Human
----	-------

技术指标

Our [Abpromise guarantee](#) covers the use of **ab90792** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

生物活性 Activity: 100 pmol of protein can bind > 80 pmol of GDP.

应用 Functional Studies

SDS-PAGE

纯度 > 90 % SDS-PAGE.

形式 Liquid

制备和贮存

稳定性和存储 Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.

pH: 8.00

Constituents: 0.00044% GDP, 0.095% Magnesium chloride, 0.0156% Beta mercaptoethanol, 0.79% Tris HCl, 0.58% Sodium chloride

This product is an active protein and may elicit a biological response in vivo, handle with caution.

常规信息

序列相似性

Belongs to the small GTPase superfamily. Rab family.

细胞定位

Cytoplasmic vesicle. Secretory granules.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors