

## Product datasheet

# Recombinant rat CNTF protein ab123755

### 概述

---

产品名称	重组大鼠CNTF蛋白
蛋白长度	Full length protein

### 描述

---

性质	Recombinant
来源	Escherichia coli

### 氨基酸序列

Accession [P20294](#)

种属 Rat

序列  
 AFAEQTPLTLHRRDLSSRSIWLARKIRSDLTALMESYVKHQGLNKNINLD  
 SVDGVPVASTDRWSEMTEAERLQENLQAYRTFQGMLTKLLEDQRVHFTPT  
 EGDHFQAIHTLMLQVSAFAYQLEELMVLLLEQKIPENEADGMPATVGDGGL  
 FEKKLWGLKVLQELSQWTVRSIHDLRVISSHQMGISALESHYGAKDKQM

分子量	23 kDa
氨基酸	2 to 200

### 技术指标

---

Our [Abpromise guarantee](#) covers the use of **ab123755** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

**生物活性** Fully biologically active by its ability to phosphorylate STAT3 in several cells lines.

**应用** SDS-PAGE

**纯度** > 95 % SDS-PAGE.  
 ab123755 is purified by proprietary chromatographic techniques. Purity by SDS-PAGE and gel filtration is greater than 99%.

**形式** Lyophilised

### 制备和贮存

---

**稳定性和存储** Shipped at 4°C. Store at -20°C.

Constituent: 0.025% Sodium carbonate

This product is an active protein and may elicit a biological response in vivo, handle with caution.

## 复溶

Centrifuge the vial prior to opening. Reconstitute in sterile ddH<sub>2</sub>O or 0.4 % NaHCO<sub>3</sub> adjusted to pH 8-9 to a concentration of 100 µg/ml. This solution can then be diluted into other aqueous buffer, preferably in the presence of carrier protein. Reconstituted CNTF should be stored at 4°C for 2-7 days and as working aliquots at -20°C for future use. Avoid freeze/thaw cycles.

## 常规信息

---

功能	CNTF is a survival factor for various neuronal cell types. Seems to prevent the degeneration of motor axons after axotomy.
组织特异性	Nervous system.
序列相似性	Belongs to the CNTF family.
细胞定位	Cytoplasm.

---

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors