

Product datasheet

Anti-RhoGDI antibody [EPR3772] ab108977

重组 RabMAb

2 References 1 图像

概述

产品名称	Anti-RhoGDI抗体[EPR3772]
描述	兔单克隆抗体[EPR3772] to RhoGDI
宿主	Rabbit
经测试应用	适用于: WB, ICC 不适用于: Flow Cyt, IHC-P or IP
种属反应性	与反应: Mouse, Rat, Human
免疫原	Synthetic peptide within Human RhoGDI aa 1-100. The exact sequence is proprietary.
阳性对照	Jurkat, HeLa, HL-60 and NIH/3T3 cell lysates
常规说明	A trial size is available to purchase for this antibody.

Our RabMAb[®] technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to [RabMAb[®] patents](#)

This product is a recombinant rabbit monoclonal antibody.

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C.
存储溶液	pH: 7.20 Preservative: 0.05% Sodium azide Constituents: 0.1% BSA, 40% Glycerol, 9.85% Tris glycine, 50% Tissue culture supernatant
纯度	Tissue culture supernatant
克隆	单克隆
克隆编号	EPR3772
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab108977** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/1000 - 1/5000. Detects a band of approximately 26 kDa (predicted molecular weight: 23 kDa).
ICC		1/50 - 1/100.

应用说明 Is unsuitable for Flow Cyt, IHC-P or IP.


靶标

功能 Regulates the GDP/GTP exchange reaction of the Rho proteins by inhibiting the dissociation of GDP from them, and the subsequent binding of GTP to them.

序列相似性 Belongs to the Rho GDI family.

细胞定位 Cytoplasm.

图片



Western blot - Anti-RhoGDI antibody [EPR3772] (ab108977)

All lanes : Anti-RhoGDI antibody [EPR3772] (ab108977) at 1/1000 dilution

Lane 1 : Jurkat cell lysates

Lane 2 : HeLa cell lysates

Lane 3 : HL-60 cell lysates

Lane 4 : NIH/3T3 cell lysates

Lysates/proteins at 10 µg per lane.

Predicted band size: 23 kDa

Observed band size: 26 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors