


Product datasheet

Anti-Secretin receptor/SCTR antibody ab85565

1 图像

概述

产品名称	Anti-Secretin receptor/SCTR抗体
描述	兔多克隆抗体to Secretin receptor/SCTR
宿主	Rabbit
经测试应用	适用于: WB, ELISA, IHC-P
种属反应性	与反应: Human 预测可用于: Mouse 
免疫原	Synthetic peptide corresponding to Human Secretin receptor/SCTR (N terminal).
阳性对照	Lung tissue.
常规说明	Protein previously labeled as SCTR.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: 0.01% Sodium azide Constituents: PBS, 50% Glycerol
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab85565** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 50 kDa.
ELISA		Use a concentration of 1 µg/ml.


应用	Ab评论	说明
IHC-P		Use a concentration of 5 µg/ml.

靶标

相关性 SCTR (Secretin receptor) is a receptor for secretin. The activity of Secretin receptor is mediated by G proteins which activate adenylyl cyclase. SCTR is an integral membrane protein and it belongs to the G-protein coupled receptor 2 family.

细胞定位 Cell Membrane; multi pass membrane protein

图片



ab85565 at 5µg/ml staining Secretin receptor/SCTR in Human lung tissue by immunohistochemistry (formalin-fixed, paraffin-embedded sections).

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Secretin receptor/SCTR antibody (ab85565)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors