


Product datasheet

Anti-SGSM3 antibody ab99422

1 图像

概述

产品名称	Anti-SGSM3抗体
描述	山羊多克隆抗体to SGSM3
宿主	Goat
经测试应用	适用于: WB
种属反应性	与反应: Human 预测可用于: Mouse, Rat, Dog 
免疫原	Synthetic peptide: RDHENYVACSRSHRR , corresponding to internal sequence amino acids 469-483 of Human SGSM3 (NP_056520.2). Run BLAST with Run BLAST with
阳性对照	A549 lysate

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: 0.02% Sodium Azide Constituents: 0.5% BSA, Tris buffered saline, pH 7.3
纯度	Immunogen affinity purified
纯化说明	ab99422 was purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab99422** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use a concentration of 2 - 3 µg/ml. Predicted molecular weight: 85 kDa.

靶标

功能	May play a cooperative role in NF2-mediated growth suppression of cells.
组织特异性	Widely expressed.
序列相似性	Belongs to the small G protein signaling modulator family. Contains 1 Rab-GAP TBC domain. Contains 1 RUN domain. Contains 1 SH3 domain.
细胞定位	Cytoplasm.

图片



Anti-SGSM3 antibody (ab99422) at 2 µg/ml +
A549 cell lysate at 35 µg

Predicted band size: 85 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors