

Product datasheet

Anti-SLIRP antibody ab51523

★★★★★ 2 Abreviews 6 References 1 图像

概述

产品名称	Anti-SLIRP抗体
描述	兔多克隆抗体to SLIRP
宿主	Rabbit
经测试应用	适用于: IHC-Fr, WB
种属反应性	与反应: Human 不与反应: Mouse
免疫原	Synthetic peptide made to a C terminal region of human SLIRP.
阳性对照	Human liver lysate

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
存储溶液	Preservative: 0.05% Sodium Azide Constituents: Tris glycine, 150mM Sodium chloride
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab51523** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
IHC-Fr	★★★★★	Use at an assay dependent concentration.


应用	Ab评论	说明
WB	★★★★★	Use a concentration of 0.5 - 2 µg/ml. Detects a band of approximately 12 kDa (predicted molecular weight: 12 kDa).

## 靶标

**相关性** SLIRP is a RNA binding protein that acts as a nuclear receptor corepressor. It may act by binding the SRA RNA, and repressing the SRA-mediated nuclear receptor coactivation. SLIRP binds the STR7 loop of SRA RNA. It is also able to repress glucocorticoid (GR), androgen (AR), thyroid (TR) and VDR-mediated transactivation. It contains 1 RRM (RNA recognition motif) domain. SLIRP is predominantly found in mitochondria with a limited role in the nucleus where it is recruited to nuclear receptor target promoters.

**细胞定位** Mitochondrial and Nuclear; Predominantly mitochondrial.

## 图片



Anti-SLIRP antibody (ab51523) at 0.6 µg/ml + human liver lysate

**Predicted band size:** 12 kDa

**Observed band size:** 12 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise,

please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors