

Product datasheet

Anti-Snf1lk2 antibody [S15G10] ab115567

1 图像

概述

产品名称	Anti-Snf1lk2抗体[S15G10]
描述	小鼠单克隆抗体[S15G10] to Snf1lk2
宿主	Mouse
经测试应用	适用于: WB, IP, IHC-P, ICC/IF
种属反应性	与反应: Human
免疫原	Recombinant Human Snf1lk2 protein (BAA34501.3).
阳性对照	Human Pancreas tissue.

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	pH: 7.20 Preservative: 0.09% Sodium azide Constituent: 99% PBS
纯度	Protein G purified
克隆	单克隆
克隆编号	S15G10
同种型	IgG1

应用

Our [Abpromise guarantee](#) covers the use of **ab115567** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明
WB		1/250. Predicted molecular weight: 104 kDa.

应用	Ab评论	说明
IP		Use at an assay dependent concentration.
IHC-P		Use a concentration of 10 µg/ml. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.
ICC/IF		1/100 - 1/400.

靶标

功能	Phosphorylates 'Ser-794' of IRS1 in insulin-stimulated adipocytes, potentially modulating the efficiency of insulin signal transduction. Inhibits CREB activity by phosphorylating and repressing TORCs, the CREB-specific coactivators.
序列相似性	Belongs to the protein kinase superfamily. CAMK Ser/Thr protein kinase family. AMPK subfamily. Contains 1 protein kinase domain. Contains 1 UBA domain.
细胞定位	Cytoplasm.

图片



ab115567, at 10µg/ml, staining Snf1k2 in Formalin-fixed, Paraffin-embedded Human Pancreas tissue by Immunohistochemistry followed by biotinylated secondary antibody, alkaline phosphatase-streptavidin and chromogen.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Promise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors