

Product datasheet

Anti-Syndecan-1 antibody [B-A38], prediluted ab1410

2 References

概述

产品名称	Anti-Syndecan-1抗体[B-A38], prediluted
描述	小鼠单克隆抗体[B-A38] to Syndecan-1, prediluted
宿主	Mouse
经测试应用	适用于: IHC-P, IHC-Fr
种属反应性	与反应: Human
免疫原	U266 cell line derived from Hodgkins disease
阳性对照	Tonsil

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: Proprietary, not sodium azide or thimerosal. Constituents: PBS, Tissue Culture Supernatant. pH 7.3; Protein carrier.
克隆	单克隆
克隆编号	B-A38
骨髓瘤	unknown
同种型	IgG1
轻链类型	unknown

应用

Our [Abpromise guarantee](#) covers the use of **ab1410** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
IHC-P		
IHC-Fr		

应用说明

IHC:

This antibody is pre-diluted for use with an ABC system.

For formalin fixed paraffin embedded sections, perform heat mediated antigen retrieval before commencing with IHC.

靶标

功能

Cell surface proteoglycan that bears both heparan sulfate and chondroitin sulfate and that links the cytoskeleton to the interstitial matrix.

序列相似性

Belongs to the syndecan proteoglycan family.

细胞定位Membrane.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors