


Product datasheet

Anti-TMEM184A antibody ab113708

2 图像

概述

产品名称	Anti-TMEM184A抗体
描述	兔多克隆抗体to TMEM184A
宿主	Rabbit
经测试应用	适用于: WB, ELISA, IHC-P
种属反应性	与反应: Rat, Human 预测可用于: Mouse 
免疫原	Synthetic peptide from the C terminus of Human TMEM184A, conjugated to KLH.
阳性对照	Rat brain tissue lysate; Human prostate tissue

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.02% Sodium azide Constituent: 99% PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab113708** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use a concentration of 1 - 2 µg/ml. Predicted molecular weight: 46 kDa.
ELISA		Use at an assay dependent concentration.


应用	Ab评论	说明
IHC-P		Use a concentration of 5 µg/ml.

## 靶标

**序列相似性** Belongs to the TMEM184 family.

**细胞定位** Membrane.

## 图片




Western blot - Anti-TMEM184 antibody (ab113708)

**All lanes :** Anti-TMEM184 antibody (ab113708) at 1 µg/ml

**Lane 1 :** Rat brain tissue lysate

**Lane 2 :** Rat brain tissue lysate with blocking peptide

**Predicted band size:** 46 kDa



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-TMEM184 antibody (ab113708)

ab113708, at 5 µg/ml, staining TMEM184 in formalin fixed, paraffin embedded Human prostate tissue by Immunohistochemistry.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

**Our Abpromise to you: Quality guaranteed and expert technical support**

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors