

## Product datasheet

# Anti-TNFRSF14 antibody [MM0332-5S1] ab89479

### 1 References

#### 概述

产品名称	Anti-TNFRSF14抗体[MM0332-5S1]
描述	小鼠单克隆抗体[MM0332-5S1] to TNFRSF14
宿主	Mouse
经测试应用	适用于: WB, IHC-P
种属反应性	与反应: Human
免疫原	Recombinant fragment corresponding to the extracellular domain of Human TNFRSF14

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Constituent: PBS
纯度	Protein G purified
纯化说明	Purified from tissue culture supernatant and 0.2 µm filter sterilised.
克隆	单克隆
克隆编号	MM0332-5S1
同种型	IgG1

#### 应用

Our [Abpromise guarantee](#) covers the use of **ab89479** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/500 - 1/2000. Predicted molecular weight: 30 kDa.
IHC-P		1/20 - 1/100.

#### 靶标

<b>功能</b>	Receptor for BTLA. Receptor for TNFSF14/LIGHT and homotrimeric TNFSF1/lymphotoxin-alpha. Involved in lymphocyte activation. Plays an important role in HSV pathogenesis because it enhanced the entry of several wild-type HSV strains of both serotypes into CHO cells, and mediated HSV entry into activated human T-cells.
<b>组织特异性</b>	Widely expressed, with the highest expression in lung, spleen and thymus.
<b>序列相似性</b>	Contains 3 TNFR-Cys repeats.
<b>翻译后修饰</b>	N-glycosylated.
<b>细胞定位</b>	Membrane.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors