

Product datasheet

TOM1L2 peptide (462-475) ab45647

概述

产品名称 TOM1L2多肽(462-475)

描述

性质 Synthetic

氨基酸序列

种属 Human

序列 C-RKKPERSEDALFAL

氨基酸 462 to 475

技术指标

Our [Abpromise guarantee](#) covers the use of **ab45647** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用 Blocking - Blocking peptide for Anti-TOM1L2 antibody ([ab24369](#))

形式 Liquid

制备和贮存

稳定性和存储 Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

常规信息

功能 Probable role in protein transport. May regulate growth factor-induced mitogenic signaling.

组织特异性 Ubiquitously expressed with higher expression in heart and skeletal muscle.

序列相似性 Belongs to the TOM1 family.

Contains 1 GAT domain.

Contains 1 VHS domain.

结构域 The GAT domain mediates interaction with TOLLIP.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors