

## Product datasheet

# Anti-TOMM70A antibody ab89624

★★★★★ 1 Abreviews 1 References 3 图像

### 概述

产品名称	Anti-TOMM70A抗体
描述	小鼠多克隆抗体to TOMM70A
宿主	Mouse
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Full length human TOMM70A protein (NP_055635).
阳性对照	Human liver tissue lysate, HepG2 cell lysate and TOMM70A transfected 293T cell lysate

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	pH: 7.20 Constituent: 2.68% PBS
纯度	Protein A purified
克隆	多克隆
同种型	IgG

### 应用

Our [Abpromise guarantee](#) covers the use of **ab89624** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 67 kDa.

### 靶标

<b>功能</b>	Receptor that accelerates the import of all mitochondrial precursor proteins.
<b>序列相似性</b>	Belongs to the Tom70 family. Contains 10 TPR repeats.
<b>细胞定位</b>	Mitochondrion outer membrane.


## 图片



Anti-TOMM70A antibody (ab89624) at 1 µg/ml  
+ human liver tissue lysate at 50 µg

**Predicted band size:** 67 kDa


**Observed band size:** 67 kDa



Anti-TOMM70A antibody (ab89624) at 1 µg/ml  
+ HepG2 cell lysate at 50 µg

**Predicted band size:** 67 kDa

**Observed band size:** 75 kDa



**All lanes :** Anti-TOMM70A antibody  
(ab89624) at 1 µg/ml

**Lane 1 :** TOMM70A transfected 293T cell  
lysate

**Lane 2 :** Non transfected 293T cell lysate

Lysates/proteins at 25 µg per lane.

**Predicted band size:** 67 kDa

**Observed band size:** 75 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors